

ORP-206A Online ORP Sensor

User Manual



YANTAI CHEMINS INSTRUMENT CO., LTD.

Tel: 0535-3463801 0571-89870583

E-mail: service@chemins-tech.com service@chemins-env.com

Website: www.chemins-tech.com

Address: No. 15, Entrepreneurship Base, Development Zone, Zhaoyuan City,

Shandong Province



User Notes

- Please read the instruction carefully before using and save it for reference.
- Please follow the instructions and precautions.
- When receiving the instrument, please open the packaging carefully, inspect equipment's damage level in case of transportation, if you found spoiled equipment, please immediately notify the manufacturer and distributor, and retain the packaging, in order to send back to processing.
- When the instrument is in trouble, please don't repair it by yourself, please directly contact the maintenance department of the manufacturer.



Content

User Notes	s2
Ι,	Application environment4
II、	Technical performance and specifications4
1.	Technical parameter4
2.	Dimensional drawing5
III、	Installation and electrical connection5
1.	Installation5
2.	Electrical connection5
IV 、	Maintenance6
1.	Use and maintenance6
2.	ORP standard solution preparation method6
3.	Calibration7
V,	Quality and service7
1.	Quality assurance7
2.	Accessories and spare parts7
3.	After-sales service commitment7
Appen	dix data communication8

$I \mathrel{\scriptstyle\diagdown}$ Application environment

Used in environmental water quality monitoring, acid/alkali/salt solution, chemical reaction process, industrial production process, it can meet the requirements of online ORP measurement for most industrial applications.

- Signal output: RS-485 (Modbus/RTU protocol).
- Convenient connection to third-party equipment such as PLC, DCS, industrial control computer, general controller, paperless recording instrument or touch screen.
- Double high-impedance differential amplifier with strong anti-interference and fast response.
- The patented ORP electrode, the internal reference solution oozes extremely slowly from the microporous salt bridge under the pressure of at least 100KPa (1Bar), and its forward bleed continues for more than 20 months. Such a reference system is very stable and the electrode life is extended by a factor of two compared to conventional industrial electrodes.
- Easy to install: 3/4 NPT pipe thread for easy submersible installation or installation in pipes and tanks.
- IP68 protection level.


$II \searrow \mbox{Technical performance and specifications}$

Model	ORP-206A		
Measuring Range	-1500 \sim +1500mV		
Resolution	1mV		
Accuracy	±6mV, ±0.5°C		
Working temperature	0∼65℃		
Working pressure	<0.1MPa		
Power supply	12~24VDC ±10%		
Signal output	RS-485(Modbus/RTU)		
Wetted material	POM		
Installation	Immersion mounting, 3/4 NPT thread		
Cable length	5 meters, other lengths can be customized		
Calibration	One point calibration		
Power consumption	0.2W@12V		
Protection grade	IP68		

1. Technical parameter




2. Dimensional drawing



Note: The sensor connector is m16-5 core waterproof connector male.

III_{\sim} Installation and electrical connection

1. Installation



Note: The sensor should not be installed upside down or horizontally when installed, at least at an angle of 15 degrees or more.

2. Electrical connection

- a) Red line power cord (12 ~ 24V)
- b) Black line ground (GND)
- c) Blue line 485A
- d) white line 485B



After wiring is completed, it should be carefully checked to avoid incorrect connections before powering up.

Cable specification: Considering that the cable is immersed in water (including sea water) for a long time or exposed to the air, all the wiring points are required to do waterproof treatment, the user cable should has certain corrosion resistance.

$I\!V\,{\scriptstyle\diagdown}\,$ Maintenance

1. Use and maintenance

When measuring the ORP electrode, it should be cleaned in distilled water (or deionized water) first, and the filter paper should be used to absorb moisture to prevent impurities from being introduced into the liquid to be tested. 1/3 of the electrode should be inserted into the solution to be tested.

The electrode should be washed when not in use, inserted into a protective cover with a 3.5 mol/L potassium chloride solution, or inserted into a container with a 3.5 mol/L potassium chloride solution.

Check if the terminal is dry. If it is stained, wipe it with absolute alcohol and dry it. Avoid long-term immersion in distilled water or protein solution and prevent contact with silicone grease. With a longer electrode, its glass film may become translucent or with deposits, which can be washed with dilute hydrochloric acid and rinsed with water. The electrode is used for a long time. When a measurement error occurs, it must be calibrated with the meter to correct it.

When the calibration and measurement cannot be performed while maintaining and maintaining the electrode in the above manner, the electrode has failed. Please replace the electrode.

2. ORP standard solution preparation method

a) ORP standard solution 256mV preparation method

2.1 g of hydrazine was weighed into 200 mL of a standard pH buffer solution having a pH of 4.00 and mixed.

b) ORP standard solution 86mV preparation method

Weigh 2.1 g of hydrazine to 200 mL of a standard pH buffer solution with a pH of 6.86 and mix.

c) ORP standard solution -40mV preparation method

Weigh 2.1 g of hydrazine to 200 mL of a standard pH buffer solution with a pH of 9.18 and mix.

Note that when preparing the pH, a total of 250 mL is prepared, and 200 mL is used for use.

pH to prepare the corresponding standard solution method:

Take 250mL of distilled water from a measuring cylinder, pour it into a beaker, add a packet of calibration powder with pH=6.86/4.00/9.18, stir evenly with a glass rod until the powder is completely dissolved, and prepare a solution with pH=6.86/4.00/9.18.



3. Calibration

Place the sensor in the prepared 86mV (or 256mV or -40mV) solution and wait for 3 to 5 minutes. After the value is stable, see if the displayed value is 86mV (or 256mV or -40mV). If not, calibration is required. See the appendix for calibration instructions.

$V \smallsetminus \ensuremath{\mathsf{Q}}$ Quality and service

1. Quality assurance

- The quality inspection department has standardized inspection procedures, advanced and perfect testing equipment and means, and strictly in accordance with the regulations, to do 72-hour aging test and stability test on the product, and not to allow one unqualified product to leave the factory.
- The receiving party directly returns the product batch with a failure rate of 2%, and all the costs incurred are borne by the supplier. The reference standard refers to the product description provided by the supplier.
- Guarantee the quantity of goods and the speed of shipment.

2. Accessories and spare parts

This product includes:

- 1 sensor
- 1 copy of the manual
- 1 certificate

3. After-sales service commitment

The company provides local after-sales service within one year from the date of sale, but does not include damage caused by improper use. If repair or adjustment is required, please return it, but the shipping cost must be conceited. Damaged on the way, the company will repair the damage of the instrument for free.



Appendix data communication

1. Data Format

The default data format for Modbus communication is: 9600, n, 8, 1 (baud rate 9600bps, 1 start bit, 8 data bits, no parity, 1 stop bit).

Parameters such as baud rate can be customized.

2. Information frame format

a) Read d	lata instruction fr	ame				
06	03	XX XX	XX	хх	XX XX	
Address	Function code	Register addres	ss Number of	f registers	CRC check code (low byte first)	
b) Read d	lata response fra	me				
06	03	xx xx>	x	ХХ	ХХ	
Address	Function code	Bytes Answei	r data	CRC check	code (low byte first)	
c) Write c	data instruction f	rame				
06	06	XX XX	XX XX	XX	XX	
Address	Function code	Register address	Write data	CRC chec	k code(low byte first)	
d) Write data response frame (same data command frame)						
06	06	XX XX	XX XX	XX	XX	
Address	Function code	Register address	Write data	CRC chec	k code (low byte first)	

3. Register address

Register address	Name	Instruction	Number of registers	Access method
		Two double-byte integers,		
40001	Measured	which are measured		Deed
(0x0000)	value	values and measured	2(4 bytes)	Read
		decimal places.		
		Calibrate in the standard		
44097	Zero	solution and write the	1 (2 hydrox)) A (with a
(0x1000)	calibration	data to the actual value of	1 (2 bytes)	Write
		the standard solution.		
44103	Zero	Returns the zero	1 (2 bytes)	Read



	(0x1006)	calibration	calibration value.		
		value			
	48195	Device address	Default address is 6, data	1 (2 bytes)	Write/
	(0x2002)	Device address	range is 1-127.		Read
			Restore calibration values		
			to factory settings, write		
	48225	factory reset	data to 0.	1 (2 bytes)	\\/rito
	(0x2020)	0x2020) factory reset	Note that the sensor must	1 (2 bytes)	Write
			be calibrated again after		
		reset before use.			

4. Command example

a) Start measurement instructions Function: Obtain the ORP value of the measuring probe, and the value of ORP is mV. Request frame: 06 03 00 00 00 02 C5 BC Response frame: 06 03 04 00 78 00 00 0C EA Example of reading:

ORP value	
00 78 00 00	

For example, the ORP value: 00 78 represents the hexadecimal reading ORP value, 00 00 means the ORP value does not have a decimal point, no temperature value, converted to a decimal value of 120.

a) Calibration instructions

The ORP sensor requires only one calibration. Before using another calibration solution, the probe must be used to clean the probe. Otherwise, it will cause a large test error.

Enter the appropriate calibration command in the prepared 86mV or 256mV standard solution.

The standard solution is 86mV: Request frame: 06 06 10 00 00 56 0C 83; Response frame: 06 06 10 00 00 56 0C 83;

The standard solution is 256mV:

Request frame: 06 06 10 00 01 00 8D 2D;

Response frame: 06 06 10 00 01 00 8D 2D.

Note: The ORP value has a negative value, such as -40mV, which is calculated as:

-{FFFF-FFD7 software test value)}=-28; decimal is -40mV.



b) Set the device ID address
Function: set the Modbus device address of the electrode;
Change the device address 06 to 01. The example is as follows
Request frame: 06 06 20 02 00 01 E3 BD
Response frame: 06 06 20 02 00 01 E3 BD

5. Error Respond

If the sensor does not correctly execute the host command, it will return the following format information:

Definition	Address	Function code	Code	CRC check
Data	ADDR	COM+80H	хх	CRC 16
Number of bytes	1	1	1	2

a) CODE: 01 –Function code error

03 – Data is wrong

b) COM: The received function code